



Email: sales@paragonmed.com
800-780-5266 - Fax 954-340-2457
Export 954-345-3990

Click for more information
<http://www.paragonmed.com/autoclave.shtml>

Elara 11 Class B Autoclave

Sterilization Cycles

The Elara 11 offers a total of 7 programs: 5 sterilization programs and 2 test programs: Bowie and Dick Test and Vacuum Test. The user can disable programs which allows for savings on validation costs.

Program	Instrument and Materials	Sterilization Temp	Drying Time	Average Cycle Time
Flash 273	Unwrapped Instruments	273°F (134°C)	1 Min.	15 Min.
Wdry 273	Wrapped, Porous and Hollow A	273°F (134°C)	20 Min.	42 Min.
NoDry 250	Wrapped Instruments	250°F (121°C)	-	35 Min.
WDry 250	Porous Wrapped and Hollow A	250°F (121°C)	20 Min.	60 Min.
Delicate 250	Porous Wrapped and Hollow A	250°F (121°C)	-	45 Min.

Technical Specifications for Elara 11

Technical Specifications	Elara11
Chamber Dimensions	11 inches (280mm) x 19.8 inches (504mm)
Chamber Volume	7.5gal (28.5L)
Number of Trays	5
Tray Dimensions (WxHxD)	6.7 inches 170mm x 0.8 inches (22mm) x 16.3 inches (442mm)
Voltage(V) Freq.(Hz)	1ph 230 V 50/60 Hz
Current(A) Power(W)	10 A 2300 W
Overall Dimensions (WxHxD)	20.8 inches (530mm) x 17.4 inches (444mm) x 25.4 inches (645mm)
Shipping Weight	184 pounds (84kg) (actual device weight is 148 pounds)
Warranty	1 year, parts and labor