



Email: sales@paragonmed.com
800-780-5266 - Fax 954-340-2457
Export 954-345-3990

Click for more information
<http://www.paragonmed.com/xray.shtml>

EVA VET Classic

Digital Veterinary Dental sensors

Specifications		
	EVA-Vet Classic Size #2	EVA-Vet Classic Size #1
Part number	9992410300	9992410400 *Discontinued
Active area	1" x 1.4" (25.8 x 36 mm)	0.8" x 1.2" (20 x 30 mm)
Capsule dimensions	31.25mm x 44.6mm x 5.25mm	25.95mm x 38.9mm x 5.25mm
Image size (unpacked)	2.1 Mb	1.3 Mb
Sensor type	Active Column MOS	
Gray shades	4096 (12 bits)	
Optical imager resolution	16.7 lp/mm	
Pixel size	30 µm x 30 µm	
Sensor cable length	2 m (~6')	
Power supply	USB powered device - no external power supply needed	
Interface	USB	
Shipping weight	Approximately 5 lbs.	
Warranty	3 Years Parts and Factory Service (US and Canada)	
Safety certifications	CSA (US and Canada), CE, IEC 60601-1-2, ISO 13485	
Recommended computer requirements*	Microsoft Windows 7 Professional, Intel Core i-7 processor, 250+GB hard drive, DVD-RW drive, video card w/ 512MB RAM, high quality monitor, USB 2.0 ports, network card, internet access for remote service.	

Minimum computer requirements*	MS Windows XP SP3 32-Bit, (Home or Professional), Pentium 4 Processor @ 1.5GHz, 2 GB RAM, USP ports, CD ROM drive, 160 GB hard drive, high quality display
--------------------------------	--

*Computer not included. While not required, Internet access will greatly improve Technical Support assistance.

It is highly recommended that you consult your IT professional for a Data Back-up system.

If budget constraints do not allow for the full recommended system, purchase the best display possible.

†For manufacturers that do not support the full 500 mA, use an external self-powered USB Hub.

Specifications are subject to change without notice.

