

## Image-Vet 4G Specifications

<b>Dimensional Data</b>	
Extension arm A = 11-13/16" (30cm)	30-7/8" (78.5 cm)
Extension arm A = 23-5/8" (60cm)*	42-3/4" (108.5 cm)
Extension arm A = 31-1/2" (80cm)	50-5/8" (128.5 cm)
<b>Total Inside Reach</b>	
Extension arm A = 11-13/16" (30cm)	56-5/16" (143 cm)
Extension arm A = 23-5/8" (60cm)*	68-1/8" (173 cm)
Extension arm A = 31-1/2" (80cm)	76" (193 cm)
<b>Parked Dimensions</b>	
Extension arm A = 11-13/16" (30cm)	26" (66 cm)
Extension arm A = 23-5/8" (60cm)*	37.5" (96 cm)
37.5" (96 cm)	45" (116 cm)
<b>Product Data</b>	
Power supply	120V ±10% 60Hz 230V ±10% 50Hz
X-ray tube	70kV, 8mA
Focal spot	0.8mm IEC 336
Focus to skin distance	20cm (7-7/8") standard
X-ray field (at collimator tip)	Diameter 60mm (2-3/8")
Duty cycle	1:32
Exposure times	0.02 to 3.2s
Shipping weight	102 Lbs (38 Kg)
Warranty	2 year limited - free phone and online support for life